Micro Arrays TestCase

| Name: TC-01 Save .cell file | | | |
|--|--------------------------------|--|--|
| Creator: Pablo Rodriguez | Creator: Pablo Rodriguez | | |
| Version: 28/04/2019 | | | |
| Description/Purpose: User can save .ce | I file in his account | | |
| Actor: Registered user | Actor: Registered user | | |
| Test Data: Saved completed [valid, error] | | | |
| Precondition: User need to log in | | | |
| Dependencies: user need a uploaded file | | | |
| Steps: | Expected Result: | | |
| 1 User need to create an account 1 User create correctly his account | | | |
| 2 User need to log in 2 User log in correctly | | | |
| 3 User need to upload a cel file | 3 user's file upload correctly | | |

| Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|--|--|
| File need to have | File name can't contain numbers (2) |
| | condition Analyzed |

| Valid prove case TC-01 Save .cell file | Classes for the condition Analyzed | Result |
|---|---------------------------------------|------------------------|
| bort.cell | (1) | Success |
| bort3.cell | (2) | Error: Name not valid. |
| | | Numbers found |

| Name: TC-01 Consulting .cell | | |
|---|--|--|
| Creator: Andres Tenesaca Burgos | | |
| Version: 28/04/2019 | | |
| Description/Purpose: user can consult | saved .cell files in DB | |
| Actor: Registered user | | |
| Test Data: Display previously saved -cel | l files | |
| Precondition: User need to log in | | |
| Dependencies: user need a uploaded file | | |
| Steps: Expected Result: | | |
| 1 User need to create an account | 1 User create correctly his account | |
| 2 User need to log in 2 User log in correctly | | |
| 3 User need to upload a cel file 3 user's file upload correctly | | |
| 4 User enters a view to display | 4 user can displat the view with all his | |
| saved .cell files | saved files | |

| Condition of entry that discusses | Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|-----------------------------------|--|--|
| Display previously saved -cell | User need to upload a file | |
| files | or previous files (1) | |

| Valid prove case TC-01 Save .cell file | Classes for the condition Analyzed | Result |
|---|---------------------------------------|---------|
| viewCell_all | (1) | Success |

| Name: TC-07 Upload .cell | | |
|---|--------------------------------|--|
| Creator: Pablo Rodriguez | | |
| Version: 28/04/2019 | | |
| Description/Purpose: user can upload .co | ell files for further analysis | |
| Actor: Registered user | | |
| Test Data: Button [working, not working] | | |
| Precondition: User need to log in | | |
| Dependencies: user need a file to upload | | |
| Steps: | Expected Result: | |
| 1. User clicks upload button and a new | 1 Successful upload | |
| window appears | | |
| 2. In the window, searchs for .cell file | | |

| Condition of entry that discusses | Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|-----------------------------------|--|--|
| Button | User upload a correct file | User upload a incorrect |
| | type (1) | file type(2) |

| Valid prove case TC-01 Save .cell file | Classes for the condition Analyzed | Result |
|---|---------------------------------------|---------|
| Button | (1) | Success |
| Bort2.png | (2) | Error |
| bort3.txt | (3) | Error |
| array.cell | (4) | Success |

| Name: TC-10 Make Query | | |
|---|--|--|
| Creator: Andrés Tenesaca | | |
| Version: 28/04/2019 | | |
| Description/Purpose: user can search for | cell with an ID in 2 external database | |
| Actor: Registered user, Unregistered use | r | |
| Test Data: DB selector ["GEOquery", "ArrayServer"], Query Id [valid, invalid] | | |
| Precondition: None | | |
| Dependencies: None | | |
| Steps: | Expected Result: | |
| 1- User has to choose which DB once to | 1- Show plot from the searched ID | |
| search | | |
| 2- then user introduces array id to | | |
| search it | | |

| Condition of entry that discusses | Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|-----------------------------------|--|--|
| DB selector | Select between two data bases(1) | |
| Query Id | User enter a correct DBID(2) | User introduce an incorrect DBID(3) |

| Valid prove case | Classes for the condition | Result |
|-----------------------|---------------------------|---------|
| TC-01 Save .cell file | Analyzed | |
| GEOquery | (1) | Success |
| GDS507 | (2) | Success |
| GSM11805AAAA | (3) | Error |

| Name: TC-11 Saving Queries | | | |
|---|---------------------------------|--|--|
| Creator: Pablo Rodriguez | | | |
| Version: 28/04/2019 | | | |
| Description/Purpose: user can store prev | rious queries in user data base | | |
| Actor: Registered user | Actor: Registered user | | |
| Test Data: Query saved completed[valid, invalid] | | | |
| Precondition: Registered user | | | |
| Dependencies: None | | | |
| Steps: Expected Result: | | | |
| - Saving successfully queried array's id 1- Query ID inserted in user data base | | | |

| Condition of entry that discusses | Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|-----------------------------------|--|--|
| Query saved completed | ID successfully save (1) | Query not introduces(2) |

| Valid prove case TC-01 Save .cell file | Classes for the condition Analyzed | Result |
|---|---------------------------------------|----------------|
| GDS507 | (1) | Success |
| No query made and save | (2) | Button enabled |
| botton available | | |

| Name: TC-12 Consulting saves queries | | |
|---|---|--|
| Creator: Andres Tenesaca | | |
| Version: 28/04/2019 | | |
| Description/Purpose: user can access pro | eviously saved queries | |
| Actor: Registered user | | |
| Test Data: Display previously saved id queries | | |
| Precondition: Registered user | | |
| Dependencies: None | | |
| Steps: Expected Result: | | |
| 1. User enters a view to display saved | 1- User can view a history of saved queries | |
| queries | | |

| Condition of entry that discusses | Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|-----------------------------------|--|--|
| Display previously saved id | Enter correct path(1) | |
| queries | | |

| Valid prove case TC-01 Save .cell file | Classes for the condition Analyzed | Result |
|---|---------------------------------------|---------|
| /viewQuery_all | (1) | Success |

| Name: TC-03 Graphics and tables | | |
|--|------------------------------------|--|
| Creator: Pablo Rodriguez | | |
| Version: 28/04/2019 | | |
| Description/Purpose: print reactive graph | nics and tools | |
| Actor: Registered user, Unregistered use | ſ | |
| Test Data: checkbox [barplot, hist, heatmap] | | |
| Precondition: None | | |
| Dependencies: None | | |
| Steps: Expected Result: | | |
| 1- User selects which type of graph or | 1- Generate a graphic of .cel data | |
| table to see | | |

| Condition of entry that discusses | Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|-----------------------------------|--|--|
| checkbox | Choose diferent types of | |
| | plots(1) | |

| Valid prove case TC-01 Save .cell file | Classes for the condition Analyzed | Result |
|---|---------------------------------------|---------|
| data.table | (1) | Success |

| Name: TC-04 Register | | |
|---|---|--|
| Creator: Andres Tenesaca | | |
| Version: 28/04/2019 | | |
| Description/Purpose: register new user | | |
| Actor: Unregistered user | | |
| Test Data: name[valid, invalid], Password | d[valid, invalid], Email [valid, invalid] | |
| Precondition: The mail introduced by th | e user can't be used by other user. | |
| Dependencies: None | | |
| Steps: | Expected Result: | |
| 1- Researcher for using the app, first | 1- Create a new user | |
| puts username | | |
| 2- then tries a password, which has to | | |
| repeat | | |
| 3- enters a unique email | | |

| Condition of entry that discusses | Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|-----------------------------------|---|--|
| Name | Max. Characters 10 digits (1) | Insert more than 10 digits (2) Insert less than 4 digits (3) Insert invalid characters (4) |
| Password | Max. Characters 10 digits and numbers (5) | No numbers inserted (6) Insert less than 8 digits (7) Repeated pasword(8) |
| Email | Max. characters 20 digits with one @ in the middle and end with dot and 2 or 3 characters more.(11) | Invalid mail(12) Incorrect format(14) Email already used(16) |

| Valid prove case TC-01 Save .cell file | Classes for the condition Analyzed | Result |
|---|---------------------------------------|--|
| Kakarot | (1) | Success |
| Esto es una frase larguísima | (2) | Error: Name not valid. Length exceeds |
| К | (3) | Error: Name not valid. Length is not min |
| @@@@ | (4) | Error: Name not valid. Invalid character. |
| Bonparapara1 | (5) | Success |
| !!!!Paaaac | (6) | Error: pass not valid. |

| | | Needs 1 number at least |
|------------------------|------|---------------------------|
| Kakarot | (7) | Error: Pass not valid. |
| | | Length is not min |
| CoNtRaSeñA1 | (8) | Error: not repeated |
| | | password |
| comolagente44@hotmail. | (11) | Success |
| com | | |
| .@yahoo.es | (12) | Error: mail not valid |
| www.yahoo.es | (14) | Error: bad email format |
| Bold_one@kamisama.com | (16) | Error: email already used |

| Name: TC-05 Login | | |
|---|----------------------|--|
| Creator: Pablo Rodriguez | | |
| Version: 28/04/2019 | | |
| Description/Purpose: Login existing user | | |
| Actor: Unregistered user | | |
| Test Data: name[valid, invalid], Password[valid, invalid] | | |
| Precondition: The user name introduced by the user can't be used by other user | | |
| Dependencies: None | | |
| Steps: | Expected Result: | |
| 1- For login, first introduce username | 1- Successful log in | |
| 2- Then the password | | |

| Condition of entry that discusses | Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|-----------------------------------|--|--|
| Name | Max. Characters 10 digits (1) Correct username (5) | Insert more than 10 digits (2) Insert less than 4 digits (3) Insert invalid characters (4) Incorrect username (6) |
| Password | Max. Characters 15 digits and numbers (9) Correct password(11) | Insert less than 8 digits (7) No numbers inserted (8) Repeated pasword(10) Incorrect password(12) Empty password(13) |

| Valid prove case | Classes for the condition | Result |
|-----------------------|---------------------------|------------------------|
| TC-01 Save .cell file | Analyzed | |
| Kakarot | (1) | Success |
| Esto es una frase | (2) | Error: Name not valid. |
| larguísima | | Length exceeds |
| K | (3) | Error: Name not valid. |
| | | Length is not min |
| @@@@ | (4) | Error: Name not valid. |
| | | Invalid character. |
| Vegetta | (5) | Success |
| Pete | (6) | Error: name not found |
| Kakarot | (7) | Error: Pass not valid. |
| | | Length is not min |
| !!!!Paaaac | (8) | Error: not repeated |
| | | password |
| Bonparapara1 | (9) | Success |
| CoNtRaSeñA1 | (10) | Error: not repeated |

| | | password |
|-------------------|------|-----------------------|
| Not_another_pass1 | (11) | Success |
| @Flyingfree4 | (12) | Success |
| [Whitespace] | (13) | Error: pass cannot be |
| | | empty |

| Name: TC-06 Change template | | |
|--|----------------------------------|--|
| Creator: Andres Tenesaca | | |
| Version: 28/04/2019 | | |
| Description/Purpose: Changing Shiny app template | | |
| Actor: registered user, Unregistered user | | |
| Test Data: Select[templates] | | |
| Precondition: None | | |
| Dependencies: None | | |
| Steps: | Expected Result: | |
| 1- User selects a template | 1- Successfully changed template | |

| Condition of entry that discusses | Classes valid for the condition Analyzed | Classes not valid for the condition analyzed |
|-----------------------------------|--|--|
| Select | Template changed successfully (1) | |

| Valid prove case | Classes for the condition | Result |
|-----------------------|---------------------------|---------|
| TC-01 Save .cell file | Analyzed | |
| Cerulean | (1) | Success |
| Cosmo | (1) | Success |
| Cyborg | (1) | Success |
| Darkly | (1) | Success |